

TREATMENT PLAN REPORT

PATIENT

Name:	LUCERO, IRMA EVARISTA	IDs:	626116, 9727
Birthdate:	Tuesday, May 25, 1948	Sex:	Femenino
Primary Oncologist:	-	Primary Oncologist ID:	-
Primary Physician:	-	Primary Physician ID:	-

PLAN

General

Plan ID:	BQT
Course ID:	C2 BQT
Course Intent:	Curativo
Plan Created:	Tuesday, May 19, 2026 10:21:33 by Valentina Morales Maldonado
Plan Last Modified:	Tuesday, May 19, 2026 12:52:01 by Valentina Morales Maldonado
EL PLAN SE HA EDITADO DESPUÉS DE QUE SE HAYA GUARDADO EN LA BASE DE DATOS.	
Technique:	INTRACAVITARY
Application Setup Type:	OTHER
Seed Collections:	5
Active sources:	5 (of 5)

Radioactive source

Treatment Time:	2.01 d (estancia sencilla máx.)
Total Air Kerma Strength:	22684.38 cGy cm²
Total Curie Seconds:	28424.56
Source Model:	varía (Aprobado)
Source Model Last Modified:	varía
Half Life:	10950.00 d
Calibration Date:	varía
Calibration Activity:	varía
Calibration Strength:	varía
Treatment Activity:	varía
Treatment Strength:	varía
Treatment Date:	Tuesday, May 19, 2026 00:00:00

Dose Prescription

Target Volume:	-
Total Dose:	2500.0 cGy (2500.0 cGy / fraction)
Number of Fractions:	1
Treatment Percentage:	100.0 %
Primary Reference Point:	NS_Applicator
Planned Total Dose:	2500.0 cGy (2500.0 cGy / fraction)
Relative Dose at Primary Reference Point:	100.0 %
Dose Calculation Medium:	Homogéneo
Dose Calculation Method:	Algoritmo: TG-43
Tiempo creación: Tuesday, May 19, 2026 12:52:12	

Dosimetric Quality Parameters

Structure	ID	Value
-----------	----	-------

Image

Image ID:	Fantoma de agua
Image Name:	Fantoma de agua BV
Image Series ID:	Fantoma de agua
Series Comment:	Fantoma de agua
Image Modality:	CT
Imaging Device:	CIMAC_Rocs
Model:	BrightSpeed S
Manufacturer:	GE MEDICAL SYSTEMS
Serial Number:	*
Patient Imaging Orientation:	Decúbito supino - céfalo caudal
Patient Treatment Orientation:	Decúbito supino - céfalo caudal
Planning Coordinate System:	IEC 61217
Image Origin:	Origen del usuario desplazamiento DICOM = (0.00cm, 0.00cm, 0.00cm)
IMAGE NOT APPROVED	

Structures

Structure Set ID:		Fantoma de agua	
ID	Type	CT Value	Approval Status
Agua	Región de control	0.0 HU	No aprobado

Seed Collections

Dwell times are displayed for a single fraction

Seed Collection ID:	M4-986	Source Model:	CDCS-M LDR
Treatment AKR:	78.39 cGy cm² / h	Calibration AKR:	163.00 cGy cm² / h
Treatment Activity:	27.28 mCi	Calibration Activity:	56.74 mCi
Treatment time:	2.01 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]		Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
2.01		0.04 4.89 1.07	78.39
Seed Collection ID:	M4-985	Source Model:	CDCS-M LDR
Treatment AKR:	73.58 cGy cm² / h	Calibration AKR:	153.00 cGy cm² / h
Treatment Activity:	25.61 mCi	Calibration Activity:	53.25 mCi
Treatment time:	2.01 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]		Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
2.01		0.10 2.94 0.51	73.58
Seed Collection ID:	J5-207	Source Model:	CDCS-J LDR
Treatment AKR:	106.76 cGy cm² / h	Calibration AKR:	191.10 cGy cm² / h
Treatment Activity:	37.16 mCi	Calibration Activity:	66.52 mCi
Treatment time:	2.01 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]		Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
2.01		0.05 0.90 0.14	106.76

Planificación de braquiterapia

Planificación de braquiterapia 15.6.8

Centro Radioterapia San Juan, San Juan

Seed Collection ID:	J5-205	Source Model:	CDCS-J LDR
Treatment AKR:	106.32 cGy cm ² / h	Calibration AKR:	190.30 cGy cm ² / h
Treatment Activity:	37.01 mCi	Calibration Activity:	66.24 mCi
Treatment time:	2.01 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal

Dwell Time [d]	Coordinates [cm] (x,y,z)	Strength [cGy cm ² / h]
2.01	-1.42 -1.09 -0.28	106.32

Seed Collection ID:	J5-204	Source Model:	CDCS-J LDR
Treatment AKR:	106.04 cGy cm ² / h	Calibration AKR:	189.80 cGy cm ² / h
Treatment Activity:	36.91 mCi	Calibration Activity:	66.06 mCi
Treatment time:	2.01 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal

Dwell Time [d]	Coordinates [cm] (x,y,z)	Strength [cGy cm ² / h]
2.01	1.52 -1.44 -0.26	106.04

Reference Points

Point ID	3D-coordinates			Fraction Dose	Total Dose	Primary Point	Volume ID
	X	Y	Z				
Ad	-2.18 cm	1.91 cm	0.35 cm	2492.0 cGy	2492.0 cGy		PTV lecho
Ai	2.25 cm	1.99 cm	0.33 cm	2500.0 cGy	2500.0 cGy		PTV lecho
Bd	-4.88 cm	1.93 cm	0.30 cm	824.6 cGy	824.6 cGy		PTV lecho
Bi	5.00 cm	1.98 cm	0.34 cm	801.5 cGy	801.5 cGy		PTV lecho
NS_Applicator	-	-	-	2500.0 cGy	2500.0 cGy	X	NS_Applicator
RECTO	-3.56 cm	-1.25 cm	-3.37 cm	838.0 cGy	838.0 cGy		Rectum
SIGMA	0.40 cm	0.51 cm	1.94 cm	3148.0 cGy	3148.0 cGy		sigma
VEJIGA	-0.95 cm	0.19 cm	1.65 cm	3266.8 cGy	3266.8 cGy		Bladder

Planificación de braquiterapia 15.6.8