

TREATMENT PLAN REPORT

PATIENT

Name:	ORREGO VERA, CARINA PATRICIA	IDs:	624903, 9685
Birthdate:	Wednesday, May 07, 1975	Sex:	Femenino
Primary Oncologist:	-	Primary Oncologist ID:	-
Primary Physician:	Gabriel Agüero	Primary Physician ID:	Agüero

PLAN

General

Plan ID:	Plan1
Course ID:	BQT
Course Intent:	
Plan Created:	Tuesday, April 21, 2026 12:25:27 by Valentina Morales Maldonado
Plan Last Modified:	Tuesday, April 28, 2026 07:25:39 by Valentina Morales Maldonado
Technique:	INTRACAVITARY
Application Setup Type:	FLETCHER_SUIT
Seed Collections:	5
Active sources:	5 (of 5)

Radioactive source

Treatment Time:	2.22 d (estancia sencilla máx.)
Total Air Kerma Strength:	20417.52 cGy cm²
Total Curie Seconds:	25584.08
Source Model:	varía (Aprobado)
Source Model Last Modified:	varía
Half Life:	10950.00 d
Calibration Date:	varía
Calibration Activity:	varía
Calibration Strength:	varía
Treatment Activity:	varía
Treatment Strength:	varía
Treatment Date:	Tuesday, April 21, 2026 00:00:00

Dose Prescription

Target Volume:	-
Total Dose:	3000.0 cGy (3000.0 cGy / fraction)
Number of Fractions:	1
Treatment Percentage:	100.0 %
Primary Reference Point:	Verificación
Planned Total Dose:	3000.0 cGy (3000.0 cGy / fraction)
Relative Dose at Primary Reference Point:	100.0 %
Dose Calculation Medium:	Homogéneo
Dose Calculation Method:	Algoritmo: TG-43
	Tiempo creación: Tuesday, April 28, 2026 07:24:25

Dosimetric Quality Parameters

Structure	ID	Value
-----------	----	-------

Image

Image ID:	Fantoma de agua
Image Name:	Fantoma de agua BV
Image Series ID:	Fantoma de agua
Series Comment:	Fantoma de agua
Image Modality:	CT
Imaging Device:	CIMAC_Rocs
Model:	BrightSpeed S
Manufacturer:	GE MEDICAL SYSTEMS
Serial Number:	*
Patient Imaging Orientation:	Decúbito supino - céfalo caudal
Patient Treatment Orientation:	Decúbito supino - céfalo caudal
Planning Coordinate System:	IEC 61217
Image Origin:	Origen del usuario desplazamiento DICOM = (0.00cm, 0.00cm, 0.00cm)
IMAGE NOT APPROVED	

Structures

Structure Set ID:	Fantoma de agua		
ID	Type	CT Value	Approval Status
Agua	Región de control	0.0 HU	No aprobado

Seed Collections

Dwell times are displayed for a single fraction

Seed Collection ID:	J3pint	Source Model:	CDCS-J LDR
Treatment AKR:	67.95 cGy cm² / h	Calibration AKR:	121.50 cGy cm² / h
Treatment Activity:	23.65 mCi	Calibration Activity:	42.29 mCi
Treatment time:	2.22 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]	2.22	Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
		1.47 -0.02 -0.29	67.95
Seed Collection ID:	J3	Source Model:	CDCS-J LDR
Treatment AKR:	64.82 cGy cm² / h	Calibration AKR:	115.90 cGy cm² / h
Treatment Activity:	22.56 mCi	Calibration Activity:	40.34 mCi
Treatment time:	2.22 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]	2.22	Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
		-2.59 -0.16 -0.35	64.82
Seed Collection ID:	J5	Source Model:	CDCS-J LDR
Treatment AKR:	106.51 cGy cm² / h	Calibration AKR:	190.30 cGy cm² / h
Treatment Activity:	37.07 mCi	Calibration Activity:	66.24 mCi
Treatment time:	2.22 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]	2.22	Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
		-0.01 -0.17 0.42	106.51

Planificación de braquiterapia

Planificación de braquiterapia 15.6.8

Centro Radioterapia San Juan, San Juan

Seed Collection ID:	J5207	Source Model:	CDCS-J LDR
Treatment AKR:	106.95 cGy cm² / h	Calibration AKR:	191.10 cGy cm² / h
Treatment Activity:	37.23 mCi	Calibration Activity:	66.52 mCi
Treatment time:	2.22 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]	2.22	Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
		-0.02 1.53 0.45	106.95
Seed Collection ID:	M2	Source Model:	CDCS-M LDR
Treatment AKR:	36.94 cGy cm² / h	Calibration AKR:	76.70 cGy cm² / h
Treatment Activity:	12.86 mCi	Calibration Activity:	26.70 mCi
Treatment time:	2.22 d (estancia sencilla máx.)	Dose calculation model:	Fuente lineal
Dwell Time [d]	2.22	Coordinates [cm] (x,y,z)	Strength [cGy cm² / h]
		0.00 3.41 0.33	36.94

Reference Points

Point ID	3D-coordinates			Fraction Dose	Total Dose	Primary Point	Volume ID
	X	Y	Z				
A	2.02 cm	2.21 cm	0.31 cm	3000.0 cGy	3000.0 cGy		PTV Utero
Ad	-2.04 cm	2.12 cm	0.32 cm	3082.9 cGy	3082.9 cGy		PTV Utero
B	4.80 cm	2.32 cm	0.33 cm	821.4 cGy	821.4 cGy		PTV Utero
Bd	-4.86 cm	2.04 cm	0.31 cm	967.9 cGy	967.9 cGy		PTV Utero
recto	-5.36 cm	-1.96 cm	-1.89 cm	654.3 cGy	654.3 cGy		Rectum
Vejiga	-1.61 cm	1.37 cm	3.08 cm	1665.7 cGy	1665.7 cGy		Bladder
Verificación	-	-	-	3000.0 cGy	3000.0 cGy	X	Fantoma

Planificación de braquiterapia 15.6.8

Final de información de geometría del sistema de imágenes

TRIAL MODE – Click here for more information